

Selective quantification of IgG directly in F12 fermentation media

Morten E. Pedersen¹, Sarah I. Gad¹, Brian Sørensen¹, Jesper Østergaard^{1,2}, and Henrik Jensen^{1,2}

¹FIDA-Tech ApS, Universitetsparken 2 (C/O University of Copenhagen), 2100 Copenhagen, Denmark,

²Department of Pharmacy, University of Copenhagen, Universitetsparken 2, 2100 Copenhagen, Denmark

Introduction

A complex sample matrix still presents a significant challenge in most biomolecular binding assays, thus requiring several purification steps for obtaining reliable assay results. This leads to a long assay development time, a need for time-consuming purification processes, and ultimately increase in operational costs. Here we present a Flow-Induced Dispersion Analysis (FIDA)-based procedure for quantifying IgG directly in F12 fermentation media, without any need for purification.

FIDA is a new capillary-based technology for quantifying biomolecular concentration directly in complex solutions (e.g. plasma) and for measuring in-solution binding under native conditions. FIDA is based on Taylor dispersion in a pressure driven flow of a ligand (termed indicator, e.g. Protein-A) interacting with the analyte of interest (e.g. an IgG antibody). The indicator alone appears small (i.e. it has a small hydrodynamic radius) when it is not bound to the antibody, but upon binding it will appear larger (i.e. the complex has a larger hydrodynamic radius). The change in apparent size forms the basis for an accurate measure of analyte concentration and interaction.

Materials & Methods

FIDalyzer instrument with 480 nm LED fluorescence detection (FIDA-Tech ApS). FIDA-coated capillary (i.d.: 75 μ m, L_T : 100 cm, L_{eff} : 88 cm). 67 mM phosphate buffer pH 7.4 was used as working buffer, Protein-A-Alexa488 as indicator (100 nM), IgG from human serum was spiked (0-

1000 nM) into 80 % F12/DMEM fermentation media and used as analyte.

Sample analysis was performed by filling the capillary with analyte, followed by injection of 39 nL indicator, which was mobilized towards the detector with analyte at 400 mbar.

Results

The apparent size (hydrodynamic radius) of Protein-A was plotted as a function of IgG concentration (see figure 1). It was observed that an increase in IgG concentration led to an increase in the apparent size of Protein-A. The unbound Protein-A had a size of 5.5 nm, but upon binding to one or two IgG molecules it increased gradually to approx. 11 nm.




Figure 1: Binding curve for the interaction between Protein-A-Alexa488 and IgG, depicted as apparent size of Protein-A-Alexa488 as function of IgG concentration, established directly in 80 % F12/DMEM fermentation media, resulting in an estimated K_d of 60 nM.

A FIDA binding isotherm was fitted to the obtained data (red line in Figure 1), and the

dissociation constant (K_d) for the protein-A-IgG complex in 80 % F12 fermentation media was estimated to be 60 nM. Furthermore, the binding curve can be used to quantify IgG in F12 fermentation media.

Conclusion

The FIDA methodology is applicable for biomolecular quantification and for measuring binding affinity directly in complex samples, thereby reducing sample purification requirements and assay development time.

For more information please contact:

FIDA-Tech ApS
Universitetsparken 2
2100 Copenhagen Ø
info@fida-tech.com
Tel.: + 45 40 30 29 93